Facilities and Resources

Microscopy, Imaging and Cytometry Resources (MICR) Core

Kamiar Moin, Ph.D., Director Jessica Back, Ph.D., Associate Director 313-577-2511 http://micr.med.wayne.edu/

Environment

Overview:

MICR is a nationally recognized state-of-the-art imaging facility with established administrative policies located in the School of Medicine at Wayne State University. With a 23 year history it is the only facility of its kind at Wayne State University. MICR offers state-of-the-art technology for the evaluation of multiple fluorophores in a single specimen utilizing confocal microscopy, flow cytometry, small animal imaging and related techniques.

Services:

Include confocal microscopy, multiphoton microscopy, conventional light microscopy; FRET and FRAP, TIRF and Atomic Force Microscopy; in vivo small animal imaging (including SPECT, CT, PET, x-ray, and optical imaging); flow cytometry and ratiometric analyses (e.g., intracellular pH and ion measurement studies), as well as three- and four- dimensional image reconstruction and quantitative measurements; cell sorting; application and study design consultation; training and workshops.

Resources:

Technical Capabilities

- Conventional light and fluorescent microscopy
- Confocal microscopy
- Multiphoton microscopy
- FRET and FRAP analysis
- High-content imaging
- Flow Cytometry
- Advanced image analysis including 3D and 4D image reconstruction and quantitative analysis
- In vivo small animal imaging including Optical, X-ray, PET, SPECT and CT
- In vivo large animal PET/CT imaging.

Instrumentation

- Zeiss LSM-410 Laser Scanning Confocal Microscope equipped with the following:
 - Zeiss Axiovert 135 inverted microscope stand
 - o 633 nm HENE laser
 - o 488/568/647 nm Multiline Argon laser
 - 364 nm Ultraviolet Argon laser
 - o Confocal imaging optics
 - Three recently upgraded photomultiplier tubes allowing simultaneous detection of three fluorochromes
 - o Fully automated programmable scanning stage
 - o Specimen incubator with temperature, humidity and CO₂ control.
 - Computer command center.
- Zeiss LSM-510 META Laser Scanning Confocal Microscope with Multiphoton capability equipped with the following:
 - Zeiss Axioplan-2 upright microscope stand
 - o 458/477/488/514 nm multiline Argon Laser.
 - o 543 nm HeNe Laser.
 - o 633 nm HeNe Laser
 - o Coherent Chameleon TiSa multiphoton near infrared laser (tunable from 720-950 nm).

- Fully automated programmable scanning stage.
- Fully automated motorized lens turret
- Two conventional fluorescent detectors, a transmitted light detector, two non-descanned detectors, and the state of the art META detector with emission finger printing allowing simultaneous detection of up to 8 fluorophores.
- Specimen incubator with temperature, humidity and CO₂ control.
- Computer command center with sophisticated imaging and ratiometric software.
- Zeiss LSM-780 Laser Scanning Confocal Microscope equipped with the following:
 - Zeiss Axioexaminer upright microscope stand with fixed stage
 - o 405 nm UV Laser
 - o 458/477/488/514 nm multiline Argon Laser.
 - o 561 nm HeNe Laser.
 - o 633 nm HeNe Laser
 - o InTune pulsed white light laser (tunable from 488-640 nm).
 - o Fully automated programmable scanning stage.
 - Two conventional fluorescent detectors, a transmitted light detector, and the cutting edge 32channel GaAsP detector with emission finger printing allowing simultaneous detection of up to 10 fluorophores.
 - Specimen incubator with temperature, humidity and CO₂ control.
 - o Computer command center with sophisticated imaging and ratiometric software.
- Zeiss LSM-800 Laser Scanning Confocal Microscope equipped with the following:
 - o Zeiss Axio Observer.Z1 microscope microscope stand with fixed stage and Definite Focus™.
 - o 405 nm UV Laser.
 - 488 nm solid state laser.
 - 561 nm solid state lase.
 - o 640 nm solid state laser.
 - o Fully automated programmable scanning stage.
 - Three GaAsP fluorescent detectors and a transmitted light detector.
 - Specimen incubator with temperature, humidity and CO₂ control.
 - o Computer command center with sophisticated imaging and ratiometric software.
- Leica TCS SP5 MP Laser Scanning Confocal Microscope with Multiphoton capability equipped with the following:
 - Leica DM6000 Advanced Inverted Microscope stand.
 - o 405 nm UV Laser.
 - 458/476/488/514 nm multiline Argon Laser.
 - o 543 nm HeNe Laser.
 - o 594 nm HeNe Laser.
 - o 633 nm HeNe Laser.
 - o Fully automated programmable scanning stage.
 - Fully automated motorized lens turret.
 - Four conventional PMT detectors, a transmitted light detector, and the two Non-descanned detectors.
 - Specimen incubator with temperature, humidity and CO₂ control.
 - o Computer command center with sophisticated imaging and ratiometric software.
- Leica TCS SP8 Laser Scanning Confocal Microscope equipped with the following:
 - Leica DMi8CEL Advanced Inverted Microscope stand.
 - o 405 nm solid-state diode Laser
 - o 488 nm solid-state diode Laser.
 - o 552 nm solid-state diode Laser.
 - o 632 nm solid-state diode Laser
 - Spectra-Physics MaiTai ps TiSa multiphoton near infrared laser (tunable from 700-1020 nm).

- Fully automated programmable scanning stage.
- o Fully automated motorized lens turret.
- o Three conventional PMT detectors, a transmitted light detector.
- o Specimen incubator with temperature, humidity and CO₂ control.
- o Computer command center with sophisticated imaging and ratiometric software.
- Zeiss Cell Observer Confocal Microscope based on spinning disk technology equipped with following:
 - Fully automated Zeiss Observer.Z1 inverted microscope stand with Definite Focus™.
 - o 488/568/633 nm solid state lasers.
 - o Yokogawa CSU-X1 microlens-based spinning disk scan-head
 - o Dual CCD cameras for simultaneous two-channel detection
 - o Fully automated programmable scanning stage
 - o Specimen incubator with temperature, humidity and CO₂ control.
 - o Computer command center with sophisticated imaging software.
- Two Zeiss ApoTome Structured Illumination Microscopes equipped with the following:
 - o Fully automated programmable scanning stage.
 - Fully automated motorized lens turret
 - Fully automated sequential 6-color detection.
 - o 3D and 4D imaging for up to 6 fluorophores.
 - o Computer command center with sophisticated imaging software.
- Zeiss Axiovert-200 inverted fluorescent photomicroscope with color SPOT CCD camera, controlled by a Power Macintosh Command Center.
- Cellomics ArrayScan V^{TI} high-content imager based on Zeiss optics with Live Module[™] and temperature control for live cell imaging.
- Bruker In-Vivo MS/FX PRO live animal multimodal optical and X-ray imager.
- Bruker In-Vivo Xtreme live animal multimodal optical, X-ray and Cherenkov radiation imager.
- Siemens Inveon SPECT/CT multimodal imager.
- Concord MicroPET R4 scanner small animal PET imager.
- GE Discovery LS4 clinical/large animal PET/CT.
- Wafergen Smart Slide System for simultaneous microscopic analysis of up to six samples under different environmental parameters. The system can accommodate any microscope.
- Cell/tissue culture hood and dual incubators.
- Amnis ImageStream Mark II imaging flow cytometer equipped with the following:
 - Six channel detection
 - o 405 nm 105mW violet excitation laser
 - 488 nm 200mW blue excitation laser
 - o 642 nm 150mW Red excitation laser
 - o 785 nm darkfield laser
 - o 60X, 40X and 20X MultiMag objective system
- Becton Dickinson LSRII equipped with four lasers operating at 488nm (blue), 405nm (violet), 561 (green) and 633nm (red). The machine has fourteen-color capability,
- Becton Dickinson FACSCanto equipped with three lasers operating at 488nm (blue), 561 (green) and 633nm (red). The FACSCanto has eight-color capability.
- SONY SY-3200 13-color Cell Sorter/Analyzer with the following:
 - o 355nm 100mW ultraviolet excitation laser

- 405 nm 60mW violet excitation laser
- 488 nm 80mW blue excitation laser
- o 561 nm 75mW blue excitation laser
- o 642 nm 70mW Red excitation laser
- Thirteen channel detection
- Computer command center
- Two SONY SH-800 Cell Sorter/Analyzers with four lasers (violet-far red) excitation for six color analysis and computer command center.
- Molecular Imaging Portal (MIP), an advanced high-speed central computer system for remote image archiving and analysis:
 - o Dell PowerEdge 2950 controller server running on Linux platform
 - 45 tb of storage capacity for image archiving
 - o Dell PowerEdge 6800 dual processor application server running Matlab and Definiens Tissue Studio™ and Developer XD™ advanced image analysis software.
 - Dell ML6000 Integrated automated tape backup system with 72 tb capacity for data backup and archiving.
- Macintosh and PC workstations with advanced image analysis software including Definiens Tissue Studio[™] and Developer XD[™], Matlab, Universal Imaging's Metamorph[™] and PerkinElmer's Volocity[™].

Scientific Expertise

Core Director: Kamiar Moin, PhD, Director, has over 30 years of experience in microscopy, imaging and cytometry thus guaranteeing stability of this vital shared resource. He received formal training in all aspects of microscopy and imaging including electron microscopy (both scanning and transmission), fluorescent and confocal microscopy. He is certified in confocal microscopy by George Washington University, Washington, D.C. Having served on several expert panels and study sections, Dr. Moin is ideally qualified to direct the research capabilities of the MICR Core and recommend scientific methodology to MICR Core users.

Associate Director: Jessica Back, PhD, has extensive experience in flow cytometry and has received cross training on several services in the Core including small animal optical imaging and microscopy. She received formal training in imaging flow cytometry and the operation of the Amnis imaging flow cytometer is one of her primary responsibilities.

Support staff are fully trained and each have in-depth experience in confocal imaging (Linda Mayernik, MS; Daniel DeSantis, BS) or flow cytometry (Eric Van Buren, BS) or in vivo imaging (Xin Lu MD; Kirk Douglas, MS). These individuals guide MICR clients in proper use of equipment, recommend techniques and resolve imaging and/or cytometry challenges.

Facility and Administration

Locations:

The MICR central imaging site is located on the medical campus and housed in approximately 1000 sq. ft of space on the 6th floor of Scott Hall. The main cytometry location is housed on the 6th floor of the Hudson Webber Cancer Research Center in approximately 1000 square feet of space. The in vivo imaging lab for small animal imaging is located in the lower level of the Elliman Building, adjacent to a central animal housing location in approximately 1500 square feet of space. Satellite locations have been established at two locations on the university main campus: for confocal imaging, located in the Biological Sciences building, and confocal imaging/cytometry, located in the Integrative Biosciences Center. Total square footage is over 4000.

6339 Scott Hall (Administrative Offices & MICR Microscopy) – 577-2511 615 HW-CRC (MICR Cytometry) – 576-8341 0125 Elliman (MICR In Vivo imaging) – 577-1207; 576-8254

2313 iBio Center (Microscopy satellite location) - 577-2511

0376 Biological Sciences (Microscopy satellite location) -- 577-2511

Access: The MICR core is available to clients from 8:30 am until 5:00 pm. After hours and weekend availability are readily provided but require prior arrangements.

Institutional Support: The MICR Core is supported by the Office of the Vice President for Research and the Karmanos Cancer Institute. In addition to subsidizing usage fees, the OVPR and KCI have and continue to contribute toward the cost of instrumentation and maintenance.

iLab Solutions: iLab Solutions is a virtual software management system that allows easy access to Wayne State University Core Services by investigators throughout the university. Cores using the iLab scheduling, invoicing and reporting system are able to identify services and associated costs and to make ordering of services available twenty-four hours a day, seven days a week. All communications between the investigator (or his/her designee) within iLab are saved, and status of completion services can be tracked. In most cases, lab results can be provided directly to the investigator online, with exporting into many formats for analysis. iLab Solutions software has been integrated with the WSU Banner financial system thus streamlining the reporting process and the expediting reimbursement for services.